

Magnetic recumbent exercise bike for home fitness. This exercise bike is ideal for lower limb muscular recovery as the tension is manually adjustable. The most important details of the workout are constantly shown on the display. The user seats in a more comfortable position compared to usual exercise bikes.



Handgrip



## TECHNICAL SPECIFICATIONS

<b>Product code</b>	<b>181030</b>
Resistance System	Magnetic
Tension Levels	8 levels manual tension
Flywheel	12 kg
Display	1 window lcd
Functions	Time / speed / distance / calories / pulse / ODO
Heart Rate Monitoring	Hand pulse
Cardio belt	Optional
Handlebar	Ergonomic, nonslip
Seat	High absorption, horizontal adjustable
Crank arms	3 pcs
Maximum User Weight	150 kg
Net Weight	49.5 kg
Gross Weight	57,6 kg
Set-Up Sizes	178 x 72.5 x 128 cm (LxWxH)
Package Dimensions	109 x 44 x 70.5 cm (LxWxH)
Transportation Wheels	Yes
Certifications	CE, ROHS

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956