

Powder coated metal structure with wooden platform with nonslip covering and convenient access ramps. Telescopic, height- and width-adjustable, anatomic handrails available in powder coated metal. The knobs prevent the loosening and movement of the handrail when used by patients.

**CODE 131120**  
**2 MT PARALLEL BAR WITH METAL HANDRAIL**



**FISIOTECH SRL**

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS

Handrail length	2 mt
Base length	2 mt
Base width	70 cm
Platform width	59 cm
Handrails height	65 – 100 cm
Distance between handrails	40 – 80 cm
Single handrail maximum load	135 kg

Handrails height adjustment system



### FISIOTECH SRL



**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## OPTIONALS

Code	Description	
131200	<p><b>PATIENT STABILISER</b></p> <p>The patient stabiliser may be used on all parallel bars with platform to support the patient during exercises. The patient stabiliser is equipped with four stoppers to be mounted on the ends of the handrails and a joining plate to connect a second patient stabilizer in case it is purchased. The double patient stabilizer helps the patient when reversing the direction.</p>	
131221	<p><b>OBSTACLES</b></p> <p>Obstacles / divider board set with magnet</p> <p>Set of eight varnished wood obstacles.</p> <p>Use it on the board to create specific rehabilitation courses.</p> <p>Dimensions: 65x13x1,5 cm</p>	

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956