

Powder coated metal structure with wooden platform with nonslip covering and convenient access ramps. Telescopic, height-adjustable, anatomic handrails available in two versions – wooden or powder coated metal. The knobs prevent the loosening and movement of the handrail when used by patients.

**CODE 131130**  
**4 MT PARALLEL BAR WITH WOODEN HANDRAIL**



**FISIOTECH SRL**

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisitech.com](mailto:info@fisitech.com) - pec : [fisitech@emailsec.it](mailto:fisitech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS

Handrail length	4 mt
Base length	4 mt
Base width	70 cm
Platform width	59 cm
Handrails height	65 – 100 cm
Distance between handrails	40 – 80 cm
Single handrail maximum load	135 kg

Handrails height adjustment system



### FISIOTECH SRL




**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisitech.com](mailto:info@fisitech.com) - pec : [fisitech@emailsec.it](mailto:fisitech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## OPTIONAL ACCESSORIES

Code	Description	
131200	<p><b>PATIENT STABILISER</b> The patient stabiliser may be used on all parallel bars with platform to support the patient during exercises. The patient stabiliser is equipped with four stoppers mountable on the ends of the handrails and a joining plate to connect a second patient stabilizer in case it is purchased. The double patient stabiliser helps the patient when reversing the direction.</p>	
131221	<p><b>OBSTACLES</b> Obstacles / divider board set with magnet Set of eight varnished wood obstacles. Use it on the board to create specific rehabilitation courses. Dimensions: 65x13x1,5 cm</p>	
131230	<p><b>CONCAVE / CONVEX PLATFORMS</b> Use them on the board to create specific rehabilitation courses. Dimensions: 290x46x8,5 cm</p>	


### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

Code	Description	
131250	<p><b>PLATFORM WITH DIVIDERS</b> Use it on the board to create specific rehabilitation courses.</p> <p><b>DIVIDERS CAN BE FILLED WITH:</b></p> <ul style="list-style-type: none"> <li>• Code 131260 Motor training sandbags dimensions: 38x20 cm</li> <li>• Code 131270 Motor training PVC bags dimensions: 38x20 cm</li> <li>• Code 131280 Wooden rollers diameter 3 cm length 20 cm</li> <li>• Code 131290 Wooden sticks diameter 1 cm length 20 cm.</li> </ul>	

**FISIOTECH SRL**

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956