



Powder coated metal structure with wooden platform with nonslip covering and convenient access ramps. Telescopic, width- and height-adjustable, anatomic, powder coated metal handrails. The knobs/ratchet handles prevent the loosening and movement of the handrail when used by patients.

CODE 1311140
6 MT PARALLEL BAR WITH STEEL HANDRAIL



FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TECHNICAL SPECIFICATIONS

Handrail length	6 mt
Base length	6 mt
Base width	70 cm
Platform width	59 cm
Handrails height	65 – 100 cm
Distance between handrails	40 – 80 cm
Single handrail maximum load	135 kg

Handrails height adjustment system



FISIOTECH SRL




Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

OPTIONAL ACCESSORIES

Code	Description	
1311200	<p>PATIENT STABILISER The patient stabiliser may be used on all parallel bars with platform to support the patient during exercises. The patient stabiliser is equipped with four stoppers mountable on the ends of the handrails and a joining plate to connect a second patient stabilizer in case it is purchased. The double patient stabiliser helps the patient when reversing the direction.</p>	
1311221	<p>OBSTACLES Obstacles / divider board set with magnet Set of eight varnished wood obstacles. Use it on the board to create specific rehabilitation courses. Dimensions: 65x13x1,5 cm</p>	
1311230	<p>CONCAVE / CONVEX PLATFORMS Use them on the board to create specific rehabilitation courses. Dimensions: 290x46x8,5 cm</p>	


FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

Code	Description	
1311250	<p>PLATFORM WITH DIVIDERS Use it on the board to create specific rehabilitation courses.</p> <p>DIVIDERS CAN BE FILLED WITH:</p> <ul style="list-style-type: none"> • Code 131260 Motor training sandbags dimensions: 38x20 cm • Code 131270 Motor training PVC bags dimensions: 38x20 cm • Code 131280 Wooden rollers diameter 3 cm length 20 cm • Code 131290 Wooden sticks diameter 1 cm length 20 cm. 	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956