


The **ADE.180** couch is a 1-section, fixed height couch. This professional couch is suitable for the Bobath method.



adjustable foot

### TECHNICAL SPECIFICATIONS

Sections	1
Under bed clearance	50 cm
Width	180 cm
Length	200 cm
Height	50 cm
Max working load, evenly distributed	200 kg
Adjustable foot	1
Weight	50 kg



### ADE COUCH MODELS

Description	Code
ADE.180	116031

#### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS – PADDING

Name of material	Polyurethane foam
Reaction to fire class	1 IM (highest)
Density	30 kg/m <sup>3</sup>
Height	40 mm



## AVAILABLE UPHOLSTERY COLOURS

When placing your order, please specify the upholstery colour.



### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: [info@fisiotech.com](mailto:info@fisiotech.com) - pec: [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS – UPHOLSTERY

### TECHNICAL DATA SHEET HAPPY AE I

<i>Resistance to penetration by body fluids (blood and urine)</i>	0-20 kPa	ISO 16603: 2004	<b>NO PENETRATION</b>
<i>Reaction ignition of upholstered furniture exposed to a small flame</i>		UNI 9175: 1987 UNI 9175/FA1: 1994	<b>CLASS 1.IM</b> <i>Item tested with flame retardant foam, density: 25-30 Kg/m<sup>3</sup>. Different foams can give different results.</i>
<i>Smouldering cigarette test -Propane flame test</i>		IMO - 2010 FTP CODE Part 8 "Test for upholstered furniture" par 7.2 - 7.3	<b>PASS</b>
<i>Specific features</i>	<i>Antibacterial                      Fungicide</i>		
<i>Contains biocidal active substances Tiabendazolo and Folpet to help preservation of the surface and oppose degradation caused by microbial growth (according to EU Regulation n. 528/2012)</i>			
<i>Kindly note: pearl colours not suitable if used under intense crackling</i>			
<i>Composition</i>	<i>PLASTICISED P.V.C.</i>		<b>88%</b>
	<i>POLYESTER</i>		<b>12%</b>

**Cleaning rules**

*Clean with damp soapy cloth and rinse well with clean water. Do not use solvents, bleaches, synthetic detergents, polishes or aerosol spray.*

**Attention**

*Normally synthetic leather light colours should not come in contact with no fixed closed clothes (such as jeans and similar) not to cause stains and rings that cannot be cleaned up.*

*The overmentioned data correspond to average values communicated in good faith, but they cannot be considered as a contract document. The product use has to be considered in terms of what has to be realized and it cannot be held up responsible when it is improperly used. Slight colour differences among several lots are admitted.*

	STANDARDS	U.M.	MINIMUM ENVIRONMENTAL CRITERIA LIMITS	OBTAINED VALUES	
<i>Phthalates</i>	UNI EN ISO 14389: 2014	%	0,1	< 0,06	
<i>Arylamines</i>	UNI EN 14362-1: 2012	mg/Kg	≤ 30	NOT DETECTABLE	
<i>Determination of formaldehyde free and hydrolysed formaldehyde</i>	UNI EN ISO 14184-1:2011	mg/Kg	≤ 75	< 3	
<i>Determination of formaldehyde free and hydrolysed formaldehyde (classroom furniture)</i>	UNI EN ISO 14184-1:2011	mg/Kg	≤ 20	< 3	
<i>Determination of metal content</i>	UNI EN 16711-2:2015	mg/Kg	ANTIMONY	≤ 30	NOT DETECTABLE
			ARSENIC	≤ 1	NOT DETECTABLE
			CADMIUM	≤ 0,1	NOT DETECTABLE
			CHROME	≤ 2	NOT DETECTABLE
			COBALT	≤ 4	NOT DETECTABLE
			COPPER	≤ 50	NOT DETECTABLE
			LEAD	≤ 1	NOT DETECTABLE
			MERCURY	≤ 0,02	NOT DETECTABLE
			NICKEL	≤ 1	NOT DETECTABLE

\* Test certificates available on request

#### FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

**GENERIC OPTIONAL ACCESSORIES**

**Code 120020 FLOOR PAPER ROLL HOLDER**

Paper roll holder for examination, massage therapy, and treatment couches.

Easy to move, powder-coated metal structure.

Dimensions: length 71 cm, width 25 cm, height 81 cm



**Code 120060 DOUBLE-PLY PAPER ROLL**

6 rolls

Dimensions: width 60 cm – length 65 m.



**FISIOTECH SRL**

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: [info@fisiotech.com](mailto:info@fisiotech.com) - pec: [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956