


ERA armchair for gynaecological examination electrically adjustable. Three sections.



CARATTERISTICHE TECNICHE

Sections	3
Width	54 cm
Length	180 cm
Height – electrically adjustable	53 - 85 cm
Max working load, evenly distributed	180 kg
Backrest angle – electrically adjustable	+90°
Footrest angle – electrically adjustable	-90°
Electric Trendelenburg position	21°
Two rigid front castors	75 mm
Two rubber stoppers	rear
Leg supports manually adjustable	in height and rotation
Stainless steel bowl	removable
Weight	90 kg



ERA COUCH MODELS

Description	Code
ERA	119034

FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TECHNICAL SPECIFICATIONS – PADDING

Name of material	Polyurethane foam
Reaction to fire class	1 IM (highest)
Density	30 kg/m ³
Height	40 mm



TECHNICAL SPECIFICATIONS – UPHOLSTERY



Origine ITALIA / Origin ITALY

Caratteristiche Properties	Norma Norm	Unita' Units	Valori Values	Tolleranza Allowance
Composizione spalmatura / Coating composition		100% PVC		
Resistenza a trazione Tensile resistance	Long. Long. Trasv. Across	UNI EN ISO 1421	N/5cm 300	minimo minimum
Allungamento a rottura Ultimate elongation	Long. Long. Trasv. Across	UNI EN ISO 1421	% 40 % 145	minimo minimum
Resistenza alla lacerazione Tear resistance	Long. Long. Trasv. Across	UNI EN ISO 4674 meth.B	N 20 N 15	minimo minimum
Solidità allo sfregamento Colour fastness to rubbing		UNI EN ISO 105-X12 : 2003 Scala dei grigi / Grey scale	Secco/Dry	4
			Umido/Wet	4
Martindale		EN 13520 : ISO 12947 EN 13520 : ISO 12947	100.000	
Classe ignifugo Flame retardant norms	UNI 9175-CLASSE 1IM UNI EN 1021-1/2	ISO 3795	BS 5852 TB 117	FMVSS 302
Formaldeide Formaldehyde	ESENTE			
Lavaggi Washability	Superficie lavabile con acqua e sapone neutro Surface washable with water and neutral soap			

UPHOLSTERY COLOURS

When placing your order, please specify the upholstery colour.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

ERA TABLE GYNECOLOGY SERIES OPTIONAL ACCESSORIES

Code 129020 PAPER ROLL HOLDER

Metal structure to be placed under the backrest of the couch.



Code 129030 HEADREST

Removable padded headrest to give the patient higher comfort.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956