



The **ERMES** couch is a four-section couch. The height can be adjusted electronically. This couch is suitable for physical therapy, massage therapy, and examination.



TECHNICAL SPECIFICATIONS

Sections	4
Under bed clearance	17 cm
Lateral offset	32 cm
Width	65 or 80 cm
Length	190 cm
Max working load, evenly distributed	200 kg
Height	47 - 80 cm
Armrests adjustment	Manual
Power supply	230 V 50/60 Hz 1,5A
Actuator thrust	8000 N
IP Code	IP 66
Headrest angle – adjustable with gas spring	+60° – 20°
Breathing hole and bung	as standard
Weight	60 kg



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

ERMES COUCH MODELS

Description	Code
Ermes with foot pedal	112031
Ermes LARGE with foot pedal	112041

HEIGHT ADJUSTMENT SYSTEMS

Ermes with 24 V foot pedal



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TECHNICAL SPECIFICATIONS – PADDING

Name of material	Polyurethane foam
Reaction to fire class	1 IM (highest)
Density	30 kg/m ³
Height	40 mm



TECHNICAL SPECIFICATIONS – UPHOLSTERY

Origine ITALIA / Origin ITALY						
Caratteristiche Properties		Norma	Norm	Unita' Units	Valori Values	Tolleranza Allowance
Composizione spalmatura / Coating composition		100% PVC				
Resistenza a trazione	Long. Long.	UNI EN ISO 1421	N/5cm	300	minimo minimum	
Tensile resistance	Trasv. Across					
Allungamento a rottura	Long. Long.	UNI EN ISO 1421	%	40	minimo minimum	
Ultimate elongation	Trasv. Across					
Resistenza alla lacerazione	Long. Long.	UNI EN ISO 4874 meth.B	N	20	minimo minimum	
Tear resistance	Trasv. Across					
Solidità allo sfregamento		UNI EN ISO 105-X12 : 2003 Scala dei grigi / Grey scale	Secco/Dry	4		
Colour fastness to rubbing			Umido/Wet	4		
Martindale		EN 13520 : ISO 12947		100.000		
		EN 13520 : ISO 12947				
Classe ignifugo		UNI 9175-CLASSE 1IM	ISO 3795	BS 5852	FMVSS 302	
Fiame retardant norms		UNI EN 1021-1/2		TB 117		
Formaldeide		ESENTE				
Formaldehyde						
Lavaggi		Superficie lavabile con acqua e sapone neutro				
Washability		Surface washable with water and neutral soap				

UPHOLSTERY COLOURS

When placing your order, please specify the upholstery colour.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

GENERIC OPTIONAL ACCESSORIES

Code 120020 FLOOR PAPER ROLL HOLDER

Paper roll holder for examination, massage therapy, and treatment couches.

Easy to move, powder-coated metal structure.

Dimensions: length 71 cm, width 25 cm, height 81 cm



Code 120040 COUCH COVER

Terry towel couch cover with elasticated hem. It is suitable for any 65 cm treatment and examination couch. Washable, blue fabric.

Code 120050 COUCH COVER

Terry towel couch cover with elasticated hem. It is suitable for any 80 cm treatment and examination couch. Washable, blue fabric.



Code 120060 DOUBLE-PLY PAPER ROLL

6 rolls

Dimensions: width 60 cm – length 65 m.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

SIMPLE SERIES OPTIONAL ACCESSORIES

Code 120010 PAPER ROLL HOLDER

Metal structure to be placed under the final section of the couch.



Code 122010 CASTORS

Four retractable castors to be used only to move the couch
Ø 75 mm.



Code 120030 COT SIDES

Cot sides ensuring the patient's safety while lying on the couch. Easy to use, with a simple manual mechanism.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956