

Fisiotek 3000 series for passive rehabilitation of lower limbs. They are intended for use in rehabilitation to restore joint movement in both surgically and nonsurgically treated medical conditions.



MODEL	CODE
Fisiotek 3000GS	134010
Fisiotek 3000TS	134020

## TECHNICAL SPECIFICATIONS

## FISIOTEK 3000GS

## FISIOTEK 3000TS

Hip and knee mobilization	✓	✓
Ankle mobilization	-	✓
Use memory card	-	✓
Speed control flexion extension	✓	✓
Working time control	✓	✓
Force control	✓	✓
Extension automatic increase	✓	✓
Flexion automatic increase	✓	✓
Pause during extension	✓	✓
Pause during flexion	✓	✓
Warm Up cycles	✓	✓
Adjustable foot support	✓	-
Knee range of movement	-10° ÷ 120°	0° ÷ 110°
Ankle range of movement	-	-20° ÷ 0° ÷ 40°
Hip range of movement	7° ÷ 115°	7° ÷ 115°
Extension automatic increase limit	✓	✓
Flexion automatic increase limit	✓	✓
Repetitions to extension limit	✓	✓
Repetitions to flexion limit	✓	✓

### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

## OPTIONAL ACCESSORIES

### Code 134030 MOBILE PROGRAMMING KEYBOARD

With its graphic, user-friendly display, the keyboard makes the equipment fully and simply programmable. The graphic symbols provided are self-explanatory.



### Code 134040 TROLLEY FOR FISIOTEK

Easy and functional, it is fitted with nonslip supports for the Fisiotek machine and with a shelf.



### Code 134050 FISIOTEK 3000 EXTENSION

This accessory can only be mounted on FISIOTEK 3000 TS. It can be used in the rehabilitation of limbs between 61 cm and 72 cm in length for a range from 0° to 110°. Using same accessory, you can reach 135° of bending in a limb from 72 to 100 cm in length.



## FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956