


Fisiotek LT is used for the passive rehabilitation of the shoulder. It is used in rehabilitation to restore joint movement in both surgically and nonsurgically treated medical conditions. This device is height-adjustable.



MODEL	CODE
Fisiotek LT	134060

TECHNICAL SPECIFICATIONS

Elevation in flexion range
Elevation in abduction
Internal-external rotation range
Speed
Range of movement

FISIOTEK LT

0° ÷ 180°
Complete physiological range
90° ÷ 0° ÷ 90°
2% sec ÷ 4.5%/sec
Adjustment by electromechanical limit switch

FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

OPTIONAL ACCESSORIES

**Code 134070 ARM SUPPORT FOR ROTATION
MOVEMENT**



Code 134080 BATTERY KIT



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956