




The **GEA** couch is a 2-section, metal-framed, fixed height couch. This couch can be used by rehabilitation and physiotherapy professionals.



## TECHNICAL SPECIFICATIONS

Sections	2
Width	65 cm
Length	190 cm
Height	78 cm
Max working load, evenly distributed	150 kg
Headrest angle – adjustable with “rastofix”	+ 60°
Breathing hole and bung	as standard
Weight	50 kg



## GEA COUCH MODELS

Description	Code
Gea fixed height	1116A10S0S

### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS – PADDING

Name of material	Polyurethane foam
Reaction to fire class	1 IM (highest)
Density	30 kg/m <sup>3</sup>
Height	40 mm



## AVAILABLE UPHOLSTERY COLOURS

 Beige P1	 Light grey G2	 Taupe P2	 Dark grey G3	 Yellow G1
 Mustard G4	 Orange A1	 Red R1	 Burgundy B2	 Violet V4
 Dark blue B3	 Blue B1	 Dark brown M1	 Earth M2	 Green V1
 Forest green V2	 Leaf green V3	 Black N1		

### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: [info@fisiotech.com](mailto:info@fisiotech.com) - pec: [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS – UPHOLSTERY

FEATURE	TEST METHOD	UNIT	LIMIT VALUE
Tensile strength	UNI EN ISO 1421 Method 1	N	Lengthwise $\geq 210$ Crosswise $\geq 90$
Elongation strength	UNI EN ISO 1421 Method 1	%	Lengthwise $\geq 85$ Crosswise $\geq 140$
Determination of tear force	EN ISO 13937-2: 2000	daN	Lengthwise $\geq 2,5$ Crosswise $\geq 2,2$
Colour fastness to crocking	UNI EN ISO 105 X12: 2003 UNI EN ISO 105 D02: 1998	grey scale	Dry 4 Bleach 4/5 Ethyl alcohol 3/4
Determination of abrasion resistance: Martindale abrader	UNI EN ISO 5470-2: 2004	cycles	Method 1 (direct) 65.000 Method 2 (reverse) 210.000
Reaction ignition Fireproof class	UNI 9175: 1987 UNI 9175/FA1: 1994		CLASSE 1.IM
Free and hydrolyzed formaldehyde	UNI EN ISO 14184-1: 2011	mg/Kg	< 3
Phthalates	UNI EN ISO 14389: 2014	%	< 0,06
Arylamines	UNI EN ISO 14362-1: 2012	mg/Kg	Not detectable
Cleaning	Clean with damp soapy cloth and rinse well with clean water. Do not use solvents, bleaches, synthetic detergents, polishes or aerosol spray.		
Specific features	Antibacterial Fingicide		
Composition	Plasticised PVC 88% - Polyester 12%		

**FIRE RESISTANCE ACCORDING TO:**  
IMO - 2010 FTP CODE Part 8  
"Test for upholstered furniture"  
par 7.2 SMOULDERING CIGARETTE TEST  
par 7.3 PROPANE FLAME TEST  
EN 1021-1: 2006 smouldering cigarette test  
EN 1021-2: 2006 match flame equivalent test

## GEA COUCH SET-CLOSE SERIES OPTIONAL ACCESSORIES

### Code 1216010 PAPER ROLL HOLDER

Metal structure to be placed under the final section of the couch.



### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: [info@fisiotech.com](mailto:info@fisiotech.com) - pec: [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## GENERIC OPTIONAL ACCESSORIES

### **Code 1210020 FLOOR PAPER ROLL HOLDER**

Paper roll holder for examination, massage therapy, and treatment couches.

Easy to move, powder-coated metal structure.

Dimensions: length 71 cm, width 25 cm, height 81 cm



### **Code 1210040 COUCH COVER**

Terry towel couch cover with elasticated hem. It is suitable for any 65 cm treatment and examination couch. Washable, grey fabric.



### **Code 1210060 DOUBLE-PLY PAPER ROLL**

6 rolls

Dimensions: width 60 cm - length 65 m.



## **FISIOTECH SRL**

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: [info@fisiotech.com](mailto:info@fisiotech.com) - pec: [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956