




The **GIOVE 5** couch is a five-section couch. The height can be adjusted electrically or hydraulically. This couch is suitable for physical therapy, massage therapy, and examination.

Thanks to its structure, height can be adjusted without lateral offset.



TECHNICAL SPECIFICATION

Sections	5
Under bed clearance	12 cm
Lateral offset	0 cm
Width	65 cm
Length	200 cm
Max working load, evenly distributed	200 kg
Height	49 – 95 cm
Power supply	230 V 50/60 Hz 1,5A
Actuator thrust	8000 N
IP Code	IP 66
Headrest angle – adjustable with gas spring	+70° –40°
Armrests adjustment	manual
Breathing hole and bung	as standard
Middle section angle	Electrical or with gas spring +20°
Footrest section angle	With gas spring +85°
Trendelenburg position	+20°
Weight	70 kg



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

GIOVE 5 COUCH MODELS

Description	Code
Giove 5 with foot pedal	118051
Giove 5 with all around foot bar switch	118052
Giove 5 hydraulic	118053
Giove 5 with aluminium all around foot bar switch	118058
Giove 5 – two motors, one used to adjust the middle section operated by hand control and one to adjust the height operated by foot pedal	118055
Giove 5 – two motors, one used to adjust the middle section operated by hand control and one to adjust the height operated by all around foot bar switch	118057
Giove 5 – two motors, one used to adjust the middle section operated by hand control and one to adjust the height operated by aluminium all around foot bar switch	118059

HEIGHT ADJUSTMENT SYSTEMS



Giove 5 with 24V foot pedal



Giove 5 with all around foot bar switch



Giove 5 with all around aluminium foot bar switch



Giove 5 hydraulic

FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

AVAILABLE STRUCTURE COLOURS



■ **BASE RAL 7016**



■ **BASE RAL 7035**

AVAILABLE UPHOLSTERY COLOURS

When placing your order, please specify the upholstery colour.



TECHNICAL SPECIFICATIONS – PADDING

Name of material	Polyurethane foam
Reaction to fire class	1 IM (highest)
Density	30 kg/m ³
Height	40 mm



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TECHNICAL SPECIFICATIONS – UPHOLSTERY

TECHNICAL DATA SHEET HAPPY AE I

<i>Resistance to penetration by body fluids (blood and urine)</i>	0-20 kPa	ISO 16603: 2004	NO PENETRATION
<i>Reaction ignition of upholstered furniture exposed to a small flame</i>		UNI 9175: 1987 UNI 9175/FA1: 1994	CLASS 1.IM <i>Item tested with flame retardant foam, density: 25-30 Kg/m³. Different foams can give different results.</i>
<i>Smouldering cigarette test -Propane flame test</i>		IMO - 2010 FTP CODE Part 8 "Test for upholstered furniture" par 7.2 - 7.3	PASS
<i>Specific features</i>	Antibacterial Fungicide		
<i>Contains biocidal active substances Tiabendazolo and Folpet to help preservation of the surface and oppose degradation caused by microbial growth (according to EU Regulation n. 528/2012)</i>			
<i>Kindly note: pearl colours not suitable if used under intense crackling</i>			
<i>Composition</i>	PLASTICISED P.V.C.		88%
	POLYESTER		12%

<i>Cleaning rules</i>	<i>Clean with damp soapy cloth and rinse well with clean water. Do not use solvents, bleaches, synthetic detergents, polishes or aerosol spray.</i>
<i>Attention</i>	<i>Normally synthetic leather light colours should not come in contact with no fixed closed clothes (such as jeans and similar) not to cause stains and rings that cannot be cleaned up.</i>

The overmentioned data correspond to average values communicated in good faith, but they cannot be considered as a contract document. The product use has to be considered in terms of what has to be realized and it cannot be held up responsible when it is improperly used. Slight colour differences among several lots are admitted.

	STANDARDS	U.M.	MINIMUM ENVIRONMENTAL CRITERIA LIMITS	OBTAINED VALUES	
<i>Phthalates</i>	UNI EN ISO 14389: 2014	%	0,1	< 0,06	
<i>Arylamines</i>	UNI EN 14362-1: 2012	mg/Kg	≤ 30	NOT DETECTABLE	
<i>Determination of formaldehyde free and hydrolysed formaldehyde</i>	UNI EN ISO 14184-1:2011	mg/Kg	≤ 75	< 3	
<i>Determination of formaldehyde free and hydrolysed formaldehyde (classroom furniture)</i>	UNI EN ISO 14184-1:2011	mg/Kg	≤ 20	< 3	
<i>Determination of metal content</i>	UNI EN 16711-2:2015	mg/Kg	ANTIMONY	≤ 30	NOT DETECTABLE
			ARSENIC	≤ 1	NOT DETECTABLE
			CADMIUM	≤ 0,1	NOT DETECTABLE
			CHROME	≤ 2	NOT DETECTABLE
			COBALT	≤ 4	NOT DETECTABLE
			COPPER	≤ 50	NOT DETECTABLE
			LEAD	≤ 1	NOT DETECTABLE
			MERCURY	≤ 0,02	NOT DETECTABLE
			NICKEL	≤ 1	NOT DETECTABLE

* Test certificates available on request

FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

GENERIC OPTIONAL ACCESSORIES

Code 120020 FLOOR PAPER ROLL HOLDER

Paper roll holder for examination, massage therapy, and treatment couches.


Easy to move, powder coated metal structure.

Dimensions: length 71 cm, width 25 cm, height 81 cm



Code 120040 COUCH COVER

Terry towel couch cover with elasticated hem. It is suitable for any 65 cm treatment and examination couch. Washable, grey fabric.



Code 120060 DOUBLE-PLY PAPER ROLL

6 rolls

Dimensions: width 60 cm – length 65 m.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TOP SERIES OPTIONAL ACCESSORIES

Code 128010 PAPER ROLL HOLDER

Metal structure to be placed under the final section of the couch.



Code 128020 CASTORS

Four retractable castors to be used only to move the couch.
Ø 75 mm.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956