




The **GIOVE 7** couch is a seven-section couch. The height can be adjusted electronically. This couch is suitable for physical therapy, massage therapy, and examination.

Thanks to its structure, height can be adjusted without lateral offset.



## TECHNICAL SPECIFICATIONS

Sections	7
Under bed clearance	12 cm
Lateral offset	0 cm
Width	65 cm
Length	200 cm
Max working load, evenly distributed	200 kg
Height	49 – 95 cm
Power supply	230 V 50/60 Hz 1,5A
Actuator thrust	8000 N
IP Code	IP 66
Headrest angle – adjustable with gas spring	+70° – 40°
Armrests adjustment	manual
Breathing hole and bung	as standard
Lateral armrests	manually foldable
Electrically operated middle section adjustment	<b>With hand control</b> +20°
Footrest adjustment	With gas spring +85°
Trendelenburg position	+ 20°
Weight	70 kg



### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

## GIOVE 7 COUCH MODELS

Description	Code
Giove 7 – two motors, one used to adjust the middle section operated by hand control and one to adjust the height operated by foot pedal	118065
Giove 7 – two motors, one used to adjust the middle section operated by hand control and one to adjust the height operated by all around foot bar switch	118067
Giove 7 – two motors, one used to adjust the middle section operated by hand control and one to adjust the height operated by aluminium all around foot bar switch	118069

## HEIGHT ADJUSTMENT SYSTEMS

Giove 7 with 24V foot pedal



Giove 7 with all around foot bar switch



Giove 7 with all around aluminium foot bar switch



### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS – PADDING

Name of material	Polyurethane foam
Reaction to fire class	1 IM (highest)
Density	30 kg/m <sup>3</sup>
Height	40 mm



## TECHNICAL SPECIFICATIONS – UPHOLSTERY



Origine ITALIA / Origin ITALY					
Caratteristiche Properties		Norma Norm	Unita' Units	Valori Values	Tolleranza Allowance
Composizione spalmatura / Coating composition		100% PVC			
Resistenza a trazione Tensile resistance	Long. Long. ----- Trasv. Across	UNI EN ISO 1421	N/5cm	300	minimo minimum
Allungamento a rottura Ultimate elongation	Long. Long. ----- Trasv. Across		%	40	minimo minimum
Resistenza alla lacerazione Tear resistance	Long. Long. ----- Trasv. Across	UNI EN ISO 4674 meth.B	N	20	minimo minimum
Solidità allo sfregamento Colour fastness to rubbing	UNI EN ISO 105-X12 : 2003 Scala dei grigi / Grey scale		Secco/Dry	4	
	EN 13520 : ISO 12947		Umido/Wet	4	
Martindale	EN 13520 : ISO 12947		100.000		
Classe ignifugo Flame retardant norms	UNI 9175-CLASSE 1IM		ISO 3795	BS 5852	FMVSS 302
Formaldeide Formaldehyde			ESENTE		
Lavaggi Washability	Superficie lavabile con acqua e sapone neutro Surface washable with water and neutral soap				

## UPHOLSTERY COLOURS

When placing your order, please specify the upholstery colour.



### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

## GENERIC OPTIONAL ACCESSORIES

### Code 120020 FLOOR PAPER ROLL HOLDER

Paper roll holder for examination, massage therapy, and treatment couches.

Easy to move, powder coated metal structure.

Dimensions: length 71 cm, width 25 cm, height 81 cm



### Code 120040 COUCH COVER

Terry towel couch cover with elasticated hem. It is suitable for any 65 cm treatment and examination couch. Washable, blue fabric.



### Code 120060 DOUBLE-PLY PAPER ROLL

6 rolls

Dimensions: width 60 cm - length 65 m.



## FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TOP SERIES OPTIONAL ACCESSORIES

### Code 128010 PAPER ROLL HOLDER

Metal structure to be placed under the final section of the couch.



### Code 128020 CASTORS

Four retractable castors to be used only to move the couch.  
Ø 75 mm.



### FISIOTECH SRL

**Operating Headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered Office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956