



Gymnic Classic used for various exercises in rehabilitation, in motor activity, for fitness and aerobics. Used in the Pilates method and in functional, sports and athletic training, too.

Medical device 93/42 EEC



### MODELS - DIMENSIONS

Diameter	Colour	Code
45 cm	Yellow	141011
55 cm	Red	141012
65 cm	Blue	141013
75 cm	Yellow	141014
85 cm	Red	141015
95 cm	Blue	141016
120 cm	Red	141017

### TECHNICAL SPECIFICATIONS

Recommended max load	120 kg
Tested max load	300 kg

#### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

**OPTIONAL ACCESSORIES**

Code	Description	
141141	Ball support "S" for balls diameter 45-55 cm	
141151	Ball support "M" for balls diameter 65-75 cm	
141160	Needle pump with flexible extension	
141161	Manual pump ideal for medium- or big-sized balls	

**FISIOTECH SRL**

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956