

Movin' Sit is an inflatable dynamic cushion that activates the trunk muscles and supports the backbone. Its shape was carefully studied in order to promote good posture by facilitating the search for a balanced pelvic position. The movement and the benefits of Movin'Sit are very similar to those of the ball. Class I Medical Device. Blue.



## MODELS – DIMENSIONS TECHNICAL SPECIFICATIONS

Length	Depth	Height	Recommended max load	Tested max load	Code
36 cm	36 cm	7 cm	120 kg	200 kg	147141

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956