



Physio Roll was born from joining two balls together and creating an equipment with a great use range and a higher stability compared to balls. It is used in the physiotherapy sector and it is ideal for those who have coordination and balance difficulties, for the disabled, for the elderly, and for nursery schools.  
Medical device 93/42 EEC



## MODELS - DIMENSIONS

Diameter	Length	Colour	Code
30 cm	50 cm	Blue	141071
40 cm	65 cm	Red	141072
55 cm	90 cm	Yellow	141073
70 cm	115 cm	Blue	141074
85 cm	130 cm	Red	141075

## TECHNICAL SPECIFICATIONS

Recommended max load	120 kg
Tested max load	300 kg

### FISIOTECH SRL



**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

**OPTIONAL ACCESSORIES**

Codice	Descrizione	
141160	Needle pump with flexible extension	
141161	Manual pump ideal for medium- or big-sized balls	

**FISIOTECH SRL**

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisitech.com](mailto:info@fisitech.com) - pec : [fisitech@emailsec.it](mailto:fisitech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956