

The **ULISSE** table is designed for examination and tilt training in supine position through a gradually increasing tilt angle. The tilt is carried out with a low voltage motor.

Equipped with four swivel castors of which 2 locking and safety straps.



TECHNICAL SPECIFICATIONS

Sections	4
Angle – electrically adjustable	0° - 90°
Lever for manual rescue	standard
Under bed clearance	15 cm
Lateral offset	0 cm
Width	65 cm
Overall width (including the lateral bars)	70 cm
Length	200 cm
Height – electrically adjusted with hand control	54 – 96 cm
Max working load	150 kg
4 castors of which 2 locking	Ø 125 mm
Safety straps	5
Independent foot boards angle – manually adjustable	+10° -23°
Power supply	230 V 50/60 Hz 1,5A
Actuator thrust	8000 N
IP Code	IP 66
Headrest angle – with gas spring	+30° -80°
Breathing hole and bung	standard
Weight	70 kg

FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY


Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

ULISSE TABLE MODELS

Description	Code
ULISSE with hand control	115064



TILT AND HEIGHT ADJUSTMENT SYSTEMS

24 V hand control for tilt and height adjustment.
The hand control can be hung to the frame.



Lever for manual rescue



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

Rev. 0 dated 05/2019

TECHNICAL SPECIFICATIONS – PADDING

Name of material	Polyurethane foam
Reaction to fire class	1 IM (highest)
Density	30 kg/m ³
Height	40 mm



TECHNICAL SPECIFICATIONS – UPHOLSTERY



Origine ITALIA / Origin ITALY					
Caratteristiche Properties		Norma Norm	Unita' Units	Valori Values	Tolleranza Allowance
Composizione spalmatura / Coating composition		100% PVC			
Resistenza a trazione Tensile resistance	Long. Long. ----- Trasv. Across	UNI EN ISO 1421	N/5cm	300	minimo minimum
Allungamento a rottura Ultimate elongation	Long. Long. ----- Trasv. Across		%	40	minimo minimum
Resistenza alla lacerazione Tear resistance	Long. Long. ----- Trasv. Across	UNI EN ISO 4674 meth.B	N	20	minimo minimum
Solidità allo sfregamento Colour fastness to rubbing			Secco/Dry	4	
				Umido/Wet	4
Martindale				100.000	
		UNI EN ISO 105-X12 : 2003 Scala dei grigi / Grey scale			
		EN 13520 : ISO 12947		EN 13520 : ISO 12947	
Classe ignifugo Flame retardant norms		UNI 9175-CLASSE 1IM	ISO 3795	BS 5852	FMVSS 302
		UNI EN 1021-1/2		TB 117	
Formaldeide Formaldehyde		ESENTE			
Lavaggi Washability		Superficie lavabile con acqua e sapone neutro Surface washable with water and neutral soap			

UPHOLSTERY COLOURS

When placing your order, please specify the upholstery colour.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

GENERIC OPTIONAL ACCESSORIES

Code 120020 FLOOR PAPER ROLL HOLDER

Paper roll holder for examination, massage therapy, and treatment couches.

Easy to move, powder-coated metal structure.

Dimensions: length 71 cm, width 25 cm, height 81 cm



Code 120060 DOUBLE-PLY PAPER ROLL

6 rolls

Dimensions: width 60 cm – length 65 m.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

ULISSE TABLE TILT SERIES OPTIONAL ACCESSORIES


Code 125060 PATIENT TABLE

The patient table can be used when the tilt table is in a supine static position. The patient table is designed so that the patient can lay its forearms on it. Its height and depth can be manually adjusted. Dimensions: 630x470x8 mm



Code 125030 AXILLARY SUPPORTS

Our axillary supports can be adjusted in height, width and depth. (set of 2 pieces)



Code 125040 PADDED PATIENT RESTRAINT

Padded band designed to provide higher and more comfortable support to the patient.



Code 125050 10 CM SAFETY STRAP

Safety strap for an easy support of the patient.
Height: 10 cm.



FISIOTECH SRL

Operating Headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered Office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956