

Chair-mounted equipment for cervical traction. Adjustable in height, fully equipped with rope, traction spreader bar, leatherette head halter and weight holder.



Product code	136014
---------------------	---------------

TECHNICAL SPECIFICATIONS

Chair legs width	50 cm
Chair legs length	65 cm
Adjustable height	150 - 225 cm

FISIOTECH SRL





Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

OPTIONAL ACCESSORIES

Code	Description	
136021	Dark brown head halter	
136041	1 kg cast iron weights. Max 8 kg per side	
136052	Ankle strap with hook	
136031	0.5 kg sandbag with eyelet	
136032	1 kg sandbag with eyelet	
136033	2 kg sandbag with eyelet	
136034	3 kg sandbag with eyelet	
136035	4 kg sandbag with eyelet	
136036	5 kg sandbag with eyelet	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956