



It is made of a hard base and a soft green dome, which is removable from the base and inflatable according to individual needs. It is a multi-functional tool, which can be used for different kinds of training, from fitness to post-traumatic knee and ankle rehab.



TECHNICAL SPECIFICATIONS

Product code	1313171
Base colour	Grey
Soft inflatable dome colour	Green
Base dimensions	Ø 47 cm
Soft inflatable dome dimensions	Ø 40 cm
Recommended max load	120 kg
Tested max load	200 kg

OPTIONAL ACCESSORIES

Code	Description	
1411160	Needle pump with flexible extension. It avoids any damage to the valve of the product. Capacity: 2 x 75 cc Length: 20,5 cm	
1411161	Manual pump to inflate medium- and big-sized balls	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956