


Powder coated metal structure with wooden platform with nonslip covering and convenient access ramps. Telescopic, height-adjustable, anatomic handrails available in two versions – wooden or powder coated metal. The knobs prevent the loosening and movement of the handrail when used by patients.

CODE 131150
CORNER PARALLEL BARS WITH STEEL HANDRAILS



FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TECHNICAL SPECIFICATIONS

Handrails length	65 – 100 cm
Platform width	59 cm
Base width	70 cm
Length of the two platforms	200-270 cm
Single handrail maximum load	135 kg

Handrails height adjustment system



FISIOTECH SRL



Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

OPTIONAL ACCESSORIES

Code	Description	
131200	<p>PATIENT STABILISER The patient stabiliser may be used on all parallel bars with platform to support the patient during exercises. The patient stabiliser is equipped with four stoppers mountable on the ends of the handrails and a joining plate to connect a second patient stabilizer in case it is purchased. The double patient stabiliser helps the patient when reversing the direction.</p>	
131221	<p>OBSTACLES Obstacles / divider board set with magnet Set of eight varnished wood obstacles. Use it on the board to create specific rehabilitation courses. Dimensions: 65x13x1,5 cm</p>	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956