

Couch stairs with chrome-plated structure and nonslip steps.



MODELS AND TECHNICAL SPECIFICATIONS

Code	Description	Step depth	Step width	First step height	Second step height	Maximum load
162101	Single step	20 cm	38 cm	22 cm	-	150 kg
162102	Two steps	20 cm	38 cm	22 cm	44 cm	150 kg

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956