

Grey varnished aluminium alloy crutches. Fixed arm support in thermoplastic material. Plastic handle with reflector. Nonslip, stain-resistant rubber stopper. 1 pair.



TECHNICAL SPECIFICATIONS CRUTCHES CODE 171010

Min/max height ground to handle	68.5-91 cm
Height handle to arm support	24 cm
Weight	0.45 kg
Max load	150 kg

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956