



Electromagnetic home fitness exercise bike. This exercise bike, also named cycle ergometer, is ideal for lower limb muscular recovery as the tension is manually adjustable. The most important details of the workout are constantly shown on the display. It is also suitable for medical stress tests.



Handgrip



TECHNICAL SPECIFICATIONS

Product code	181020
Resistance System	Electromagnetic power 8V
Tension Levels	16 levels electronically tension
Flywheel	10 kg
Display	1 colour window LCD
Functions	Time / speed / distance / calories / RPM / program diagram / effort / watt
Programs	Manual – 10 pre set (16 levels)
Heart Rate Monitoring	Hand pulse + Wireless receiver
Heart Rate Chest Belt	Optional
Handlebar	Ergonomic, multiple handles, adjustable
Seat	Graduated compression, horizontally and vertically adjustable
Maximum User Weight	120 kg
Net Weight	28 kg
Gross Weight	33 kg
Set-Up Sizes	103 x 52 x 135 cm (LxWxH)
Package Dimensions	86 x 32 x 64.5 cm (LxWxH)
Transportation Wheels	Yes
Certifications	CE, ROHS, EN957 CE

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956