

Wall mounted white varnished steel sheet. It has cold fluorescent lamps (6500 K) with an unbreakable opaline, screen-printed, acrylic diffuser.



MODELS AND TECHNICAL SPECIFICATIONS

Code	Dimensions	Distance
152030	68x31x9 cm	3 mt

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956