

Equipment for cervical traction.  
Adjustable in height, fully equipped with rope, traction spreader bar, leatherette head halter and weight holder.

Traction unit not supplied



## MODELS - DIMENSIONS

Models	Base width	Base length	Adjustable height	Code
With weight holder	62 cm	62 cm	150 - 225 cm	136012
For traction unit	62 cm	62 cm	150 - 225 cm	136013

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

## OPTIONAL ACCESSORIES

Code	Description	
136021	Dark brown head halter	
136041	1 kg cast iron weights. Max 8 kg per side	
136052	Ankle strap with hook	
136031	0.5 kg sandbag with eyelet	
136032	1 kg sandbag with eyelet	
136033	2 kg sandbag with eyelet	
136034	3 kg sandbag with eyelet	
136035	4 kg sandbag with eyelet	
136036	5 kg sandbag with eyelet	

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956