

Blue PVC-covered compressed foam mattress. 600 gr/m³, antibacterial, anti-dust mite, fire resistant. Removable upholstery, with zip and handles for easy transportation. Suitable for gym exercises.




Agglomerated foam mattress with anti-mould antibacterial fire-resistant covering. With non-slip underlayer. Removable upholstery, with zip and handles for easy transportation. Suitable for gym exercises.



MODELS – DIMENSIONS TECHNICAL SPECIFICATIONS

Code	Length	Width	Height	Density	Colour
144011	200 cm	100 cm	5 cm	25 kg / m ³	Blue
144012	200 cm	100 cm	10 cm	21 kg / m ³	Blue
144013	200 cm	100 cm	5 cm	80 kg / m ³	Blue

OPTIONAL ACCESSORIES FOR MATTRESSES

Code	Description	
144901	Horizontal mat trolley. BLUE varnished steel structure with four rubber swivel castors diameter 100 mm.	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956