

Powder coated metal structure. Telescopic, height-adjustable, anatomic handrails available in two versions – wooden or powder coated metal. The knobs prevent the loosening and movement of the handrail when used by patients.

Code 131010



Code 131020



FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TECHNICAL SPECIFICATIONS

Handrails length	3 mt
Handrails height	69 – 100 cm
Distance between handrails	60 cm
Maximum width	70 cm
Single handrail maximum load	135 kg

Handrails height adjustment system



FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956