


Exclusive polyurethane integral skin foam formulation (without HFC and CFC) with ODP (Ozone Depletion Potential) blowing agents. The completely open cell structure gives the mat excellent shock absorption features, allowing the normal airflow from every single layer of the product, the homogeneity in thickness and density in any spot of the mat. The special skin on the top and on the bottom lends further mechanical resistance. These mats guarantee great resistance to abrasions and an excellent grip.



MODELS – DIMENSIONS TECHNICAL SPECIFICATIONS

Code	Length	Width	Height	Density	Colour
144031	200 cm	100 cm	3 cm	PK 14	Blue
144032	200 cm	100 cm	4 cm	PK 16	Blue
144033	200 cm	100 cm	5 cm	PK 20	Blue

OPTIONAL ACCESSORIES FOR MATTRESSES

Code	Description	
144901	Horizontal mat trolley. BLUE varnished steel structure with four rubber swivel castors diameter 100 mm.	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956