

Wall bar with 13 oval clear-varnished beech rungs. Powder-coated steel uprights with bars to fix the product to the wall. The wall bar is useful in carrying out joint and muscular exercises, in particular preventive and therapeutic exercises against backache.

Other uses of the wall bar are rehabilitation and work out exercises, joint stretching, and exercises to correct bad posture and joint movement exercises in children.



## TECHNICAL SPECIFICATIONS

Height	Width	Code
240 cm	90 cm	142010

### FISIOTECH SRL



**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## OPTIONAL ACCESSORIES

Code	Description	
142500	Black curved backrest. Padded, upholstered in black, washable, and fire-resistant leatherette. Dimensions: Height 65 cm Width 40 cm - Depth 18 cm	
142510	Black wall bar board. Padded board upholstered in black, washable, and fire-resistant leatherette. With feet-holding belt. Dimensions: Height 180 cm – Width 30 cm	

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956