




4mm-thick mirrors with accident prevention adhesive film. In case the mirror breaks, the fragments will not fall.

Mirrors with safety film are mounted on a polystyrene board and a white MDF board.

An aluminium frame binds together the components.

The mirrors can either be mounted on a wall or on a moving inclination-adjustable trolley.

Rectangles dimensions: 10x5 cm



MODELS AND TECHNICAL SPECIFICATIONS

Code	Model	Length	Width	Thickness	Weight
143010	With grid	200 cm	100 cm	4 mm	30 kg
143020	Plain	200 cm	100 cm	4 mm	30 kg
143030	With grid	170 cm	100 cm	4 mm	30 kg
143040	Plain	170 cm	100 cm	4 mm	30 kg

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

OPTIONAL ACCESSORIES

Code	Description	
143500	Moving, inclination-adjustable trolley for 200 cm mirrors.	
143510	Moving, inclination-adjustable trolley for 170 cm mirrors.	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956