

**Stainless steel 18/10 multi-purpose trolley with two shelves.  
With two pushing handles, four swivel castors.  
Available with or without drawers.**



## MODELS AND TECHNICAL SPECIFICATIONS

Models	Length	Depth	Height	Diameter castors	Code
Without drawers	82 cm	52 cm	82 cm	125 mm	164101
With two horizontal drawers	82 cm	52 cm	82 cm	125 mm	164111

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956