

With steel frame for height and inclination adjustment. Wall-mountable with wall plugs. Steel uprights and dried beech rungs, 30 mm diameter, of which three protrude laterally. Plywood, 8mm-thick central backrest board. Width: 48cm with rubber stoppers. Height: 3 m.



TECHNICAL SPECIFICATIONS

Description	Code
Straight orthopaedic ladder	142100
Curved orthopaedic ladder	142110

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956