


Device for the plantar analysis of feet.
Varnished beech structure, crystal footboard and double LED lighting on two sides.



TECHNICAL SPECIFICATIONS CODE 151010

Type of analysis	Visual
Footboard dimensions	44x39 cm
Dimensions	50x48x35 h cm
Maximum load	135 kg
Lighting	LED
Thickness of the tempered crystal	12 mm
Mirror with polished edges	
LED absorption	2 X 3 W
Insulating class	1



Power supply 230V 50Hz
I1 Omnipolar switch
F1 and F2 fuse 6.3 x 32, 500 mA
L1 and L1 LED absorption 2 x 3 W

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956