

Balance training board made of a wooden table with non-slip covering and a hemispherical ball with adjustable inflation. The two parts are linked by hook and loop fasteners. It is effective for rehabilitation but its range of use comprises training and ankle strengthening, too.



TECHNICAL SPECIFICATIONS

Code	133090
Dimensions	Ø 37
Height	18 cm
Recommended max load	120 kg

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956