

Varnished wood rocker boards with nonslip covering for balance training and rehabilitation.



TECHNICAL SPECIFICATIONS

Max working load	120 kg
------------------	--------

MODEL DIMENSIONS

Product code	Model	Length	Width	Height
1313010	With half-circle	170 cm	60 cm	14 cm
1313020	With half-circle	80 cm	60 cm	13 cm
1313030	With half-circle	60 cm	40 cm	10 cm
1313040	With half-circle	40 cm	35 cm	9 cm
1313050	Central fulcrum	60 cm	40 cm	10 cm

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956