

The pedal exerciser is suitable for those with motor difficulties. It is suitable to recover upper and lower limb coordination and muscle tone, too.

Colour: BLACK



TECHNICAL SPECIFICATIONS



Frame	Powder coated steel
Colour	Black
Non-slip feet	4
Pedals	With straps
Friction	Adjustable with belt
Base length	47 cm
Base depth	25.5 cm
Max. height	35 cm

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956