

Pedal exerciser to recover coordination and upper and lower limb muscle tone. Four non-slip tips, pedals with adjustable foot straps, handle for effort adjustment. With display.



TECHNICAL SPECIFICATIONS

Product Code 1314520



Frame	powder coated tubular steel frame
External dimensions	49x37x34 cm
Minimum height (pedals in horizontal position)	27 cm
Maximum height (pedals in vertical position)	34 cm
Weight of the product	2,3 kg

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956