

Wooden balance training board with nonslip covering.
The cylinder moves **safely** thanks to the protruding
wooden parts on the extremities of the board that
prevent the cylinder from rolling out.



TECHNICAL SPECIFICATIONS

Product code	133080
Length	60 cm
Width	30 cm
Recommended max load	120 kg

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956