

Professional instrument with clock-like dial. Aluminium and painted metal structure. Non-toxic, nonslip rubber platform. With nonslip feet.



MODELS AND TECHNICAL SPECIFICATIONS

| Code | Description | Platform dimensions | Dial dimensions | Division | Max load |
|--------|----------------|---------------------|-----------------|----------|----------|
| 153010 | With column | 27x40x8.5 cm | Diameter 18 cm | 1000 g | 150 kg |
| 153020 | Without column | 27x40x8.5 cm | Diameter 18 cm | 1000 g | 150 kg |

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956