



Shoulder and elbow rehabilitation instrument and for shoulder and elbow mobility maintenance exercises.

The instrument is made of a metal structure with a base and a wheel, which is adjustable in height to conform to the patient's needs and to the requirement of the exercise.

Exercises are carried out either sitting or standing.

The therapist can set the difficulty of the exercises by adjusting the resistance of the wheel.



MODELS DIMENSIONS – TECHNICAL SPECIFICATIONS

Code	135050
Height	160 cm
Width	110 cm
Depth	66 cm
Base dimensions	90x60 cm

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956