

Single section clear-varnished beech wall bar with oval rungs. One of the two upper rungs is placed inwards while the other outwards to carry out suspension exercises. With zinc-plated steel joints to fix the wall bar to the wall. The wall bar is useful in carrying out joint and muscular exercises, in particular preventive and therapeutic exercises against backache. Other uses of the wall bar are rehabilitation and workout, joint stretching, and exercises to correct bad posture and joint movement exercises in children.



TECHNICAL SPECIFICATIONS

Height	Width	Code
250 cm	90 cm	142030

FISIOTECH SRL







Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

OPTIONALS

Code	Description	
142500	Black curved support Padded, upholstered in washable fire-resistant leatherette. Dimensions: Height 65 cm Width 40 cm - Depth 18 cm	
142510	Black wall bar board. Wall bar board. Padded board upholstered in black, washable, and fire-resistant leatherette. With feet-holding belt. Dimensions: Height 180 cm – Width 30 cm	
142520	Suspension frame. Suspension frame for wall bar, powder-coated steel with grid for accessories. Useful for cervical traction and pulley therapy. Dimensions: Height 40 cm – Width 64 cm Depth 87 cm.	
137114	Pulley	
136052	Ankle strap with hook	
137112	Rope length 20 m	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956