



Powder coated steel structure, zinc-chromate sliding parts. Telescopic, height-adjustable, powder coated metal anatomic handrails. The knobs prevent the loosening and movement of the handrail when used by patients. Nonslip PVC covered closed steps in two different heights. The stairs can be set in a straight position as well as in a corner position.

**CODE 1312030**  
**DOUBLE FLIGHT REHABILITATION STAIRS**



**FISIOTECH SRL**

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY


**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS

Handrails adjustment	60 – 100 cm
Step width	56 cm
Step depth	30 cm
Step height	13 – 16 cm
Single handrail maximum load	135 kg



### FISIOTECH SRL


**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## OPTIONAL ACCESSORIES

Code	Description	
1312200	Inclined ramps dimensions: 170x49 cm	

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisitech.com](mailto:info@fisitech.com) - pec : [fisitech@emailsec.it](mailto:fisitech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956