

Powder coated steel structure, zinc-chromate sliding parts. Telescopic, height-adjustable, powder coated metal anatomic handrails. The knobs prevent the loosening and movement of the handrail when used by patients. Nonslip PVC covered closed steps in two different heights. The stairs can be set in a straight position as well as in a corner position.

CODE 132030
DOUBLE FLIGHT REHABILITATION STAIRS



FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY


Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

TECHNICAL SPECIFICATIONS

Handrails adjustment	60 – 100 cm
Step width	56 cm
Step depth	30 cm
Step height	13 – 16 cm
Single handrail maximum load	135 kg



FISIOTECH SRL


Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisitech.com - pec : fisitech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956

OPTIONAL ACCESSORIES

Code	Description	
132200	Inclined ramps dimensions: 170x49 cm	

FISIOTECH SRL

Operating headquarters: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Registered office: Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisitech.com - pec : fisitech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956