

**Foldable aluminium alloy walker with four rubber stoppers.**



**Foldable aluminium alloy walker with two rubber stoppers and two rigid castors.**



## TECHNICAL SPECIFICATIONS

With four stoppers	Code 172010
With two stoppers and two castors	Code 172020
Underarm height adjustment	80 - 95 cm
Length	45 cm
Width	55 cm
Weight	2.5 kg
Max load	110 kg

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956