



Shoulder and elbow rehabilitation instrument and for shoulder and elbow mobility maintenance exercises.

The instrument is made of a wall-mounted metal structure with wheel, which is adjustable in height to conform to the patient's needs and to the requirement of the exercise.

Exercises are carried out either sitting or standing.

The therapist can set the difficulty of the exercises by adjusting the resistance of the wheel.



**MODELS  
DIMENSIONS – TECHNICAL SPECIFICATIONS**

Code	135040
Height	160 cm
Width	110 cm
Depth	30 cm
Weight	29 kg
Wall plugs	<b>NOT SUPPLIED</b>

**FISIOTECH SRL**

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : info@fisiotech.com - pec : fisiotech@emailsec.it

VAT Number IT 03765100262 - Rea TV 296956