


Wall mounted white varnished steel sheet instrument with stainless steel film holder. It has cold fluorescent lamps (6500° K) with unbreakable opaline Plexiglass diffuser.



## MODELS AND TECHNICAL SPECIFICATIONS

Code	Dimensions
152010	90x60 cm
152020	70x42 cm

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956